

PCTWORLD INTELLECTUAL PROPERTY ORGANIZATION
International Bureau

INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

(51) International Patent Classification ⁵: C12Q 1/18 // G01N 33/04, 33/12, (C12Q 1/18, C12R 1:07) (C12Q 1/18, C12R 1:46)	A1	(11) International Publication Number: WO 94/18343 (43) International Publication Date: 18 August 1994 (18.08.94)
(21) International Application Number: PCT/EP94/00359 (22) International Filing Date: 9 February 1994 (09.02.94) (30) Priority Data: 93200387.4 11 February 1993 (11.02.93) EP (34) Countries for which the regional or international application was filed: NL et al. (71) Applicant (for all designated States except US): GIST-BROCADES N.V. [NL/NL]; Wateringseweg 1, P.O. Box 1, NL-2600 MA Delft (NL). (72) Inventors; and (75) Inventors/Applicants (for US only): VAN RIJN, Ferdinand, Theodorus, Jozef [NL/NL]; Valkenlaan 7, NL-2623 HM Delft (NL). BEUKERS, Robert [NL/NL]; Sijtwinde 159, NL-2631 GZ Nootdorp (NL). KERKHOF, Johannes, H., P., M. [NL/NL]; Westereerf 2, NL-3181 JD Rozenburg (NL). (74) Agents: MATULEWICZ, Emil, Rudolf, Antonius et al.; Gist-Brocades N.V., Patents and Trademarks Department, Wateringseweg 1, P.O. Box 1, NL-2600 MA Delft (NL).		(81) Designated States: AU, BB, BG, BR, BY, CA, CN, CZ, FI, HU, JP, KP, KR, KZ, LK, LV, MG, MN, MW, NO, NZ, PL, RO, RU, SD, SK, UA, US, UZ, VN, European patent (AT, BE, CH, DE, DK, ES, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE), OAPI patent (BF, BJ, CF, CG, CI, CM, GA, GN, ML, MR, NE, SN, TD, TG). Published <i>With international search report.</i> <i>Before the expiration of the time limit for amending the claims and to be republished in the event of the receipt of amendments.</i>
(54) Title: UNIT FOR THE DETECTION OF RESIDUES OF ANTIBACTERIAL COMPOUNDS IN LIQUIDS		
(57) Abstract		
<p>The detection of residues of antibacterials such as antibiotics and sulpha compounds in liquids such as milk, water, meat juice, serum or urine is disclosed. A test unit comprises an agar medium inoculated with a suitable test organism and two or more redox indicators.</p>		
<p>Foreign Cite No. 7 Appl No. 10/578935 August 9, 2006 Robert S. Salter 0656-032US3A</p>		

FOR THE PURPOSES OF INFORMATION ONLY

Codes used to identify States party to the PCT on the front pages of pamphlets publishing international applications under the PCT.

AT	Austria	GB	United Kingdom	MR	Mauritania
AU	Australia	GE	Georgia	MW	Malawi
BB	Barbados	GN	Guinea	NE	Niger
BE	Belgium	GR	Greece	NL	Netherlands
BF	Burkina Faso	HU	Hungary	NO	Norway
BG	Bulgaria	IE	Ireland	NZ	New Zealand
BJ	Benin	IT	Italy	PL	Poland
BR	Brazil	JP	Japan	PT	Portugal
BY	Belarus	KE	Kenya	RO	Romania
CA	Canada	KG	Kyrgyzstan	RU	Russian Federation
CF	Central African Republic	KP	Democratic People's Republic of Korea	SD	Sudan
CG	Congo	KR	Republic of Korea	SE	Sweden
CH	Switzerland	KZ	Kazakhstan	SI	Slovenia
CI	Côte d'Ivoire	LI	Liechtenstein	SK	Slovakia
CM	Cameroon	LK	Sri Lanka	SN	Senegal
CN	China	LU	Luxembourg	TD	Chad
CS	Czechoslovakia	LV	Latvia	TG	Togo
CZ	Czech Republic	MC	Monaco	TJ	Tajikistan
DE	Germany	MD	Republic of Moldova	TT	Trinidad and Tobago
DK	Denmark	MG	Madagascar	UA	Ukraine
ES	Spain	ML	Mali	US	United States of America
FI	Finland	MN	Mongolia	UZ	Uzbekistan
FR	France			VN	Viet Nam
GA	Gabon				

UNIT FOR THE DETECTION OF RESIDUES OF ANTIBACTERIAL
COMPOUNDS IN LIQUIDS

The present invention relates to a method for the detection of residues of antibacterial compounds in liquids. The invention also relates to a unit for the detection of residues of antibacterial compounds in liquids and the use
5 of the unit.

Similar tests have been described in GB-A-1467439, EP-A-0005891 and DE 3613794. These documents all deal with "ready to use" tests that make use of a test organism and will give a result generally between 2½ to 3½ hours by the
10 change of colour of an acid-base or redox indicator added to the test system. The principle is that when an antibacterial compound is present in the sample liquid in a sufficient concentration to inhibit the growth of the test organism the colour of the indicator will stay the same, whilst when no
15 inhibition occurs, the growth of the test organism is accompanied by the formation of acid or reduced metabolites that will change the colour of the indicator. In all these tests a single indicator is used to detect the formation of acid or reduced metabolites.

20 According to this invention a marked reduction in test duration, up to one hour, can be achieved when a combination of two or more indicators is used. Such a shortened test duration is of importance to the user because the quality of the sample liquid is known more quickly thus
25 allowing an earlier delivery or processing etc.

Therefore the invention provides a method for detecting antibacterials in a test sample which comprises

- 2 -

a) bringing a test organism and at least two redox indicators into an agar medium,

b) allowing the test sample to come into contact with the agar medium such that antibacterials in the test sample inhibit the test organism in the agar medium.

The present invention also provides a unit for detecting antibacterials which comprises an agar medium comprising a test organism, optionally a separate nutrient source and two or more redox indicators which each can be present in the agar medium, in the test sample or in the separate nutrient source.

The unit of the present invention is useful for detecting residues of antibacterials such as sulpha compounds and antibiotics.

The unit may be used for detecting antibacterials in liquids, for example milk, water, meat juices, serum or urine.

The test organism is preferably a strain of Bacillus or Streptococcus. A preferred species of Streptococcus is Streptococcus thermophilus, more preferably Streptococcus thermophilus T101 (DSM 4022, deposited on March 3, 1987). A strain of this species may be introduced into the agar medium, preferably in concentrations of 10^5 to 10^9 colony forming units (CFU) per ml agar medium.

A preferred species of Bacillus is Bacillus stearothermophilus, more preferably Bacillus stearothermophilus var. calidolactis C953. Bacillus stearothermophilus may be introduced into the agar medium preferably in concentrations of 10^5 to 10^9 CFU per ml of agar medium.

Examples of units useful for the purpose of the invention are transparent tubes, single or in a set or combined to a block of translucent material provided with a number of holes shaped therein.

A large variety of redox indicators may be used according to the process of the present invention. Such redox indicators are also known as redox mediators, redox

- 3 -

catalysts and electron carriers. Examples of such compounds are Brilliant Black, Methylene Blue, Toluidine Blue, Safranin O, Indigo Carmin, Thionin, Gallocyanine, Nile Blue A, Brilliant Crocein MOO, Acid Yellow 38, Acid Orange 51, 5 Acid Blue 120, Basic Blue 3, Azure A, Azure B, Congo Red, 1-10 Phenanthroline, Janus Green B, Brilliant Cresyl Blue. Other redox indicators (redox mediators, redox catalysts and electron carriers) may be used as well. Such compounds are commercially available see e.g. 'Stains, Dyes and 10 Indicators', Catalogue of Aldrich Chemie. Preferably one of the indicators should give a colour change in the visible part of the spectrum. Preferred combinations are

- a) Brilliant Black and Methylene Blue
- b) Brilliant Black and Toluidine Blue and
- 15 c) Brilliant Black and Nile Blue A.

Nutrients are added to enable the multiplication of the test organism.

The unit of the present invention optionally comprises at least part of the nutrients which are not 20 incorporated in the agar medium and thus are added as a separate source, for example as a tablet or paper disc, which may be placed on the agar medium before carrying out the test. Nutrients may be present both in the agar medium and as a separate source. At least one of the redox 25 indicators may be included in the separate nutrient source.

Nutrients and one or both redox indicators, e.g. in a tablet, may also be included in the units beforehand whereby measures are preferably taken to avoid moisture transport from the agar medium into the tablet. This may be 30 done, e.g. by coating the tablet with a moisture-resistant layer, for example a wax, which coating must degrade or melt during the test procedure. A wax having a melting temperature of 35 to 55°C, preferably 40 to 45°C, is suitable.

Strain C953 of Bacillus stearothermophilus var. 35 calidolactis has been deposited with the Laboratory of Microbiology of the Technical University of Delft under the accession number LMD 74.1 in 1974 and with the Centraal

- 4 -

Bureau voor Schimmelcultures (CBS), Baarn under the accession number CBS 760.83 in 1983 where the strain is available to the public. This microorganism is very sensitive to penicillins and other antibiotics and is a fast growing microorganism. It has the additional advantage that the optimal growing temperature is high (between 50-70°C). Only a few microbial species are able to grow at this temperature. There is therefore little possibility that organisms present in the test liquid or which have otherwise been included in the test system would affect the test result.

When the test organism is a Bacillus strain, it is preferably incorporated into the agar medium in the form of a spore suspension which may be prepared according to known methods (GB-A-1467439). The spore suspension is incorporated into the agar medium by known methods (GB-A-1467439).

According to a preferred embodiment of the present application the sensitivity of the test organism is adjustable. The sensitivity may be altered by various means, for example by adding certain substances, by changing the test conditions such as the pH or concentration of buffering substances, agar or spores or by varying ratio of the volumes in the volumes of agar and test liquid. Examples of substances that may be added are nucleosides, such as adenosine, or antifolates, such as trimethoprim, ormethoprim and tetroxoprim, which improve the sensitivity of the test organism to sulpha compounds, salts of oxalic acid or hydrofluoric acid which improve the sensitivity to tetracyclines, and cysteine to diminish the sensitivity to penicillins.

It is preferred to carry out the process of the present invention in such a way that the test organism stays alive but does not multiply in the agar medium. This is generally achieved by depriving the organism of nutrients until the test is carried out or/and by maintaining the culture at a sufficiently low temperature, for example below 30°C.

- 5 -

In the detection of residues of antibacterial compounds in fluids, preferably biological fluids, such as milk, water, meat juice, serum and urine, using the units as defined herein, a predetermined amount of the sample to be tested, for example 0.05 ml to 0.5 ml is placed on the agar medium (for example 0.2 ml to 3 ml), and the contents of the unit are incubated at or near the optimal temperature for the test organism for example 63°C to 66°C during a predetermined period, for example 60 to 120 minutes, after which the indicator colour is observed, indicating the presence or absence of antibacterials above a certain minimum concentration. The test is very simple to carry out, so that the person that performs the test does not have to be specially qualified. The test is completed in 1 to 2 hours after starting the test, which is markedly shorter than other microbial test systems where only one indicator is used.

All patent applications and patents mentioned in this application are herein incorporated by reference to the same extent as if each individual application or patent was specifically and individually indicated to be incorporated by reference.

The embodiments of the present application are illustrated by means of the following examples.

Example 1

Preparation of test tubes to detect antibiotics

A solution was made of 12 g agar and 9 g sodium chloride in 1000 ml distilled water. To this solution 2.5 ml of a 0.09 M triethanolamine buffer solution (pH 8.0) was added. The final solution was sterilized for 20 minutes at 121°C and cooled to about 60°C. To this sterile solution a sufficient amount of a suspension of Bacillus stearothermophilus var. calidolactis spores in distilled water was added to give a final concentration between 10^9 and 10^{10} spores per litre and an amount of a sterile solution of

- 6 -

Brilliant Black to give a final concentration of 80 mg per litre. Sterile tubes with a diameter of about 9 mm were filled with 0.3 ml of the agar solution under aseptic conditions and immediately sealed e.g. with an aluminium foil. The contents of the test tubes was allowed to solidify while the tubes were held in an upright position. The thus prepared tubes were stored at a temperature between 5°C and 15°C.

10 Example II

Preparation of a test tube to detect antibiotics and sulpha compounds

The procedure described in Example I was followed except that together with the buffer solution an amount of a trimethoprim solution was added to give a final concentration of 60 µg per litre.

Example III

20 Preparation of nutrient tablets

A mixture was made of 100 g dextrose, 160 g Avicel PH101, 50 g tryptose, 10 g phytone peptone, 15 g precirrol and 500 mg of Toluidine Blue dissolved in 50 ml of ethyl alcohol. This mixture was sufficient to prepare 30000 tablets with a diameter of 3 mm and a thickness of 1.2 mm.

Example IV

Carrying out a test

30

A test tube, prepared according to Example I or II, was opened by puncturing the seal and a nutrient tablet prepared according to Example III was added. Of the sample, e.g. a milk sample, to be investigated, 0.1 ml was added to the test tube and the test tube was placed in an incubator (waterbath or block heater) kept at 64°C. Observations were made after 1 hour and 20 minutes to 1 hour and 40 minutes.

- 7 -

If at this time the colour of the agar medium is yellow, the sample does not contain a detectable amount of an antibacterial compound (e.g. 0.003 I.U. or less of penicillin G or 0.1 μ g or less of sulfamethazine per ml).

5 If, however, the colour of the agar medium is blue, the sample contains at least a detectable amount of an antibacterial compound (e.g. 0.006 I.U. or more of penicillin G or 0.2 μ g or more of sulfamethazine per ml).

10 An in-between concentration, thus representing the just detectable amount may give a colour of the agar medium between yellow and blue.

Example V

Comparison of the two indicator test unit with a single
15 indicator test tube

Test tubes were prepared according to Example I or II. Nutrient tablets were prepared according to Example III. Similar nutrient tablets were prepared but without Toluidine
20 Blue.

Tests were carried out according to Example IV with both types of nutrient tablets and an antibiotic-free milk sample. The test tubes in combination with the nutrients that do not contain Toluidine Blue took about one hour more
25 to change colour, that is two hours 20 minutes to 2 hours 40 minutes, when compared with the combination test tube + nutrient tablet containing Toluidine Blue.

Example VI

30 Preparation of variations and carrying out tests therewith

Test tubes were made according to Example I with the distinction that instead of a solution of Brilliant Black a solution was used containing a similar concentration of one
35 of the following redox indicators: Brilliant Crocein MOO, Acid Yellow 38, Acid Yellow 51, Acid Blue 120 or Congo Red. Nutrient tablets were made according to Example II with the

- 8 -

distinction that instead of Toluidine Blue a similar amount was used of one of the following redox indicators: Safranine O, Indigo Carmine, Thionin, Nile Blue A, Azure A, Azure B, Janus Green B, Brilliant Cresyl Blue ALD or Methylene Blue.

- 5 The test was carried out according to Example IV with test tubes and nutrient tablets prepared according to the description given above.

- A resulting change in colour indicates that the milk sample investigated does not contain detectable amounts of
10 antibiotic and/or sulpha compounds. If the colour of the test does not change the milk sample does contain such residues. The test duration varied with the chosen combination of redox indicators but was markedly shorter than that of a one indicator test, in a similar way as
15 described in Example V.

- 9 -

Claims

1. A method for detecting antibacterials in a test sample which comprises

5 a) bringing a test organism and at least two redox indicators into an agar medium,

b) allowing the test sample to come into contact with the agar medium such that antibacterials in the test sample inhibit the test organism in the agar medium.

10

2. A method according to claim 1 wherein the test organism is a strain of Bacillus or Streptococcus.

3. A method according to claim 1 or 2 where the
15 test organism is a strain of Bacillus stearothermophilus var. calidolactis or a strain of Streptococcus thermophilus.

4. A method according to any one of the preceding claims wherein the agar medium is buffered.

20

5. A method according to any one of the preceding claims wherein a nutrient source is added to the agar medium.

25 6. A method according to any one of the preceding claims wherein the nutrient source or part thereof is added in the form of a separate nutrient source such as a tablet or a paper disc.

30 7. A method according to claim 6 wherein the separate nutrient source further comprises at least one redox indicator, preferably at least two redox indicators.

8. A method according to any one of the preceding
35 claims wherein at least one of the redox indicators gives a colour change in the visible part of the spectrum.

- 10 -

9. A method according to any one of the preceding claims wherein the sensitivity of the test organism to anti bacterial compounds is adjustable.

5 10. A method according to any one of the preceding claims which further comprises adding to the agar medium at least one substance which changes the sensitivity of the test organism to antibacterial compounds.

10 11. A method according to claim 10 wherein the substance which changes the sensitivity is an antifolate.

 12. A method according to claim 11 wherein the substance which changes the sensitivity is trimethoprim,
15 ormethoprim or tetroxoprim.

 13. A method according to any one of the preceding claims wherein 10^5 to 10^9 colony forming units of test organism are present per millilitre agar medium.

20 14. A unit for detecting antibacterials which comprises an agar medium comprising a test organism, optionally a separate nutrient source and two or more redox indicators which each can be present in the agar medium or
25 the separate nutrient source.

 15. A unit according to claim 14 wherein the test organism is a strain of Bacillus or Streptococcus.

30 16. A unit according to claim 15 wherein the test organism is Bacillus stearothermophilus or Streptococcus thermophilus.

 17. A unit according to any one of claims 14 to 16
35 wherein the agar medium is buffered.

- 11 -

18. A unit according to any one of claims 14 to 17 wherein the agar medium further comprises at least one nutrient.

5 19. A unit according to claim 14 wherein the separate nutrient source comprises at least one redox indicator, preferably at least two redox indicators.

20. A unit according to any one of claims 14 to 19
10 wherein at least one of the redox indicators gives a colour change in the visible part of the spectrum.

21. A unit according to any one of claims 14, 19 or 20 wherein the separate nutrient source is in the form of a
15 tablet or paper disc.

22. A unit according to any one of claims 14 to 21 wherein the sensitivity of the test organism is adjustable.

20 23. A unit according to any one of claims 14 to 21 which further comprises at least one substance which changes the sensitivity of the test organism to antibacterial compounds.

25 24. A unit according to claim 23 wherein the substance which changes the sensitivity of the test organism is an antifolate.

25 25. A unit according to claim 24 wherein the
30 substance which changes the sensitivity of the test organism is trimethoprim or tetroxoprim.

26. A unit as claimed in any one of claims 14 to 25 when used in the method as claimed in any one of claims 1 to
35 13.

INTERNATIONAL SEARCH REPORT

International Application No

PCT/EP 94/00359

A. CLASSIFICATION OF SUBJECT MATTER

IPC 5 C12Q1/18 //G01N33/04,G01N33/12,(C12Q1/18,C12R1:07),(C12Q1/18,
C12R1:46)

According to International Patent Classification (IPC) or to both national classification and IPC

B. FIELDS SEARCHED

Minimum documentation searched (classification system followed by classification symbols)

IPC 5 C12Q

Documentation searched other than minimum documentation to the extent that such documents are included in the fields searched

Electronic data base consulted during the international search (name of data base and, where practical, search terms used)

C. DOCUMENTS CONSIDERED TO BE RELEVANT

Category*	Citation of document, with indication, where appropriate, of the relevant passages	Relevant to claim No.
Y	EP,A,0 005 891 (GIST-BROCADES N.V.) 12 December 1979 cited in the application see the whole document ---	1-26
Y	EP,A,0 322 591 (ABBOTT LABORATORIES) 5 July 1989 see page 3, line 1 - line 22 see page 4, line 16 - line 17 ---	1-26
A	DE,A,36 13 794 (MÜLLER) 29 October 1987 cited in the application see the whole document ---	1,8-12
A	EP,A,0 285 792 (VALIO MEIJERIEN KESKUSOSUUSLIIKE) 12 October 1988 see abstract -----	2,3

☐ Further documents are listed in the continuation of box C.☒ Patent family members are listed in annex.

* Special categories of cited documents:

- "A" document defining the general state of the art which is not considered to be of particular relevance
- "E" earlier document but published on or after the international filing date
- "L" document which may throw doubts on priority claim(s) or which is cited to establish the publication date of another citation or other special reason (as specified)
- "O" document referring to an oral disclosure, use, exhibition or other means
- "P" document published prior to the international filing date but later than the priority date claimed

- "T" later document published after the international filing date or priority date and not in conflict with the application but cited to understand the principle or theory underlying the invention
- "X" document of particular relevance; the claimed invention cannot be considered novel or cannot be considered to involve an inventive step when the document is taken alone
- "Y" document of particular relevance; the claimed invention cannot be considered to involve an inventive step when the document is combined with one or more other such documents, such combination being obvious to a person skilled in the art.
- "&" document member of the same patent family

Date of the actual completion of the international search

20 May 1994

Date of mailing of the international search report

14. 06. 94

Name and mailing address of the ISA

European Patent Office, P.B. 5818 Patentaan 2
NL - 2280 HV Rijswijk
Tel. (+ 31-70) 340-2040, Tx. 31 651 epo nl,
Fax: (+ 31-70) 340-3016

Authorized officer

Ceder; 0

INTERNATIONAL SEARCH REPORT

Information on patent family members

International Application No

PCT/EP 94/00359

Patent document cited in search report	Publication date	Patent family member(s)	Publication date
EP-A-0005891	12-12-79	NL-A- 7806086 AU-B- 528996 AU-A- 4774079 CA-A- 1136031 JP-A- 54159295	07-12-79 19-05-83 13-12-79 23-11-82 15-12-79
EP-A-0322591	05-07-89	AU-A- 2584588 JP-A- 1187097	01-06-89 26-07-89
DE-A-3613794	29-10-87	NL-A- 8700746	16-11-87
EP-A-0285792	12-10-88	AU-B- 609794 AU-A- 1415888 CA-A- 1316442 DE-A- 3871216 JP-A- 63269999 US-A- 4929546	09-05-91 06-10-88 20-04-93 25-06-92 08-11-88 29-05-90

THIS PAGE BLANK (USPTO)